

**FOR YOUR PROTECTION, DO NOT USE IF INNER TAMPER EVIDENT SEAL IS BROKEN OR MISSING.**

**Keep out of reach of children.**

**Preparation:** Add enough bottled water to cover powder and quickly swirl vigorously to prevent sticking of powder to the bottle. Then add bottled water to the **top of the label** and shake well for 20 seconds. Let sit **on its side** about an hour and reshake – repeat as necessary (or let sit overnight) until all is in solution. Reshake and use.

**Suggested use:** as a dietary supplement to support joint health. Use product within 90 days from preparation.

**Dosage for Companion Animals:**

1/3 teaspoon small dogs & cats  
2/3 teaspoon medium dogs  
2 teaspoons large dogs  
Equine 2 Tablespoons or 1 oz. per day

To maximize potency, product should be stored closed at room temperature and away from direct sunlight. May be refrigerated.

Date Prepared: \_\_\_\_\_

**\*\*These statements have not been evaluated by US FDA. This product is not intended to diagnose, treat, cure, or prevent any disease.**

# Maximum Wellness For Companion Animals



## Hyaluronan for Life Support

DIETARY SUPPLEMENT

**15 Fluid Ounces (444 mL)  
(when reconstituted)**

**PATENT FORMULA 10849923**

*Brought to you by our team of  
physicians, pharmacists and chemical engineers.*  
[www.HA-TOPSHELF.com](http://www.HA-TOPSHELF.com)

### SUPPLEMENT FACTS

Serving Size: 2 teaspoons (10 mL)

Servings Per Container: 45 doses

Ingredients	Amount (per Serving)	% DV
Sodium Hyaluronan (sodium hyaluronate)	5 mg 100 mg	<2 *
Proprietary Blend	200 mg	*
D-mannitol, sorbitol, cyclodextrin, colloidal silver.		

\*No Daily Value has been established.

Other ingredients: none.

**\*\*Our neutral TopShelf HA is a high molecular weight patented, Hyaluronic acid product that provides the highest daily dosage for joint and health maintenance. The powder-to-liquid formula guarantees the freshest, most stable HA on the market.**

Manufactured for:  
HA-TOP SHELF®  
1330 Commercial Street  
Warsaw, MO 65355 USA



Please recycle

