

FOR YOUR PROTECTION, DO NOT USE IF INNER TAMPER EVIDENT SEAL IS BROKEN OR MISSING.

Keep out of reach of children.

If you are under a doctor's care or taking medication, you should consult with your doctor before taking any supplement.

Preparation: Add enough bottled water to cover powder and quickly swirl vigorously to prevent sticking of powder to the bottle. Then add bottled water **one inch above the top of the label** and shake well for 20 seconds. Let sit **on its side** about an hour and reshake -- repeat as necessary (or let sit overnight) until all is in solution. Reshake and use.

Suggested use: Take two teaspoons per day (humans) as a dietary supplement to support joint health. First time users should double the dosage for the first 7-10 days to replenish lost HA. Use product within 90 days from preparation.

Dosage for Humans:

2 teaspoons (150 doses)

To maximize potency, product should be stored closed at room temperature and away from direct sunlight. May be refrigerated.

Date Prepared: _____

****These statements have not been evaluated by US FDA. This product is not intended to diagnose, treat, cure, or prevent any disease.**

Brought to you by our team of physicians, pharmacists and chemical engineers.

Maximum Wellness For the Human Race



DIETARY SUPPLEMENT
PATENT FORMULA 10849923

HYALURONAN for joint and life support

Powder-to-Liquid Gel Formulation
(for optimum HA stability & effectiveness)

JUST ADD BOTTLED WATER

50 Fluid Ounces (1.5 Liters)
When Reconstituted

SUPPLEMENT FACTS

Serving Size: 2 teaspoons (10 mL)
Servings Per Container: 150 doses

Ingredients	Amount (per Serving)	% DV
Sodium Hyaluronan (sodium hyaluronate)	5 mg 100 mg	<2 *
Proprietary Blend	200 mg	*
D-mannitol, sorbitol, cyclodextrin, colloidal silver.		
*No Daily Value has been established.		
Other ingredients: none.		

****Our neutral TopShelf HA is a high molecular weight patented, Hyaluronic acid product that provides the highest daily dosage for joint and health maintenance. The powder- to-liquid formula guarantees the freshest, most stable HA on the market.**



Please recycle



Manufactured for:
HA-TOP SHELF®
1330 Commercial Street
Warsaw, MO 65355 USA
www.ha-topshelf.com



PL-DK1365-01-50